

#### About Samsung Electronics Co. Ltd.

 $Samsung\ Electronics\ Co.\ Ltd.\ is\ a\ global\ leader\ in\ technology, opening\ new\ possibilities\ for\ people\ everywhere.\ Through\ relentless\ innovation$  $and\ discovery, we are\ transforming\ the\ worlds\ of\ televisions,\ smartphones,\ personal\ computers,\ printers,\ cameras,\ home\ appliances,\ LTE\ systems,$ medical devices, semiconductors and LED solutions. We employ 236,000 people across 79 countries with annual sales of US\$187.8 billion. To discover more, please visit www.samsung.com



For more information, please write to hme.marketing@samsung.com Corporate office: 2nd, 3rd & 4th Floor, Tower-C, Vipul Tech Square, Sector – 43, Golf Course Road, Gurgaon-122002, Haryana (India). Visit: www.samsung.com/in/healthcare

# **SAMSUNG ULTRASOUND H60**

### PERFORMANCE IN STYLE







#### **PERFORMANCE IN STYLE**

The new H60 delivers revolutionary performance while adhering to the fundamental medical equipment design principles of simplicity, lightness and versatility that make it easy to operate in compact spaces. Features like a 10.1-inch touch screen, 18.5-inch LED monitor and slim body personify Samsung's devotion to practical and usable ideas. The H60 is also equipped with a cutting-edge hybrid imaging engine and S-Vue Transducer technology to provide high-resolution images.

#### **Hybrid Imaging Engine**

The enhanced hardware and newly added software engine of the H60 produces the clearest and most accurate images through innovative and intelligent technology.




#### **S-Vue Transducer**





The extremely responsive S-Vue Transducer provides broader bandwidth and higher sensitivity. This results in superb image resolution even in technically challenging cases. In addition, the

ergonomically designed and lightweight Transducer ensures users experience less fatigue. The CV1-8AD is approximately 30% lighter than the conventional model.



# **UPHOLDING THE HIGHEST ACCURACY STANDARDS**

High image quality is the key to an effective ultrasound procedure. The H60 does just that and beyond, utilising innovative technologies to deliver an unparalleled performance in high-level imaging.


#### 18.5-inch Wide LED Monitor

The H60's 18.5-inch LED ultrasound monitor provides high-quality images with wide viewing angles and the additional benefit of low power consumption.



#### ClearVision™

The noise reduction filter improves edge enhancement and creates sharper 2D images for optimal diagnostic performance. The integration of Samsung's specialised technology results in a notable improvement of image quality.




Early foetus with ClearVision™

#### S-Flow™

S-Flow<sup>™</sup> helps in detecting peripheral blood vessels through advanced colour flow imaging with superior sensitivity. The H60's S-Flow<sup>™</sup> facilitates scanning even in technically difficult diagnostic situations leading to increased accuracy and resulting in better treatment.



Umbilical cord with S-Flow™



#### **CRYSTAL CLEAR IMAGING**


Expect nothing less than sharp and high quality images, thanks to the innovative technologies powering our diagnostic systems. Clear images and powerful performance are the signature qualities of our ultrasound machines.



Foetal heart with S-Flow™



BPD and FL measurement



Foetal abdomen ascites with Multi-Slice View™



Gestational sac with e-Motion Marker™



Foetal abdomen circulation





Foetal abdomen ascites in Oblique View™




Early foetus



Umbilical cord in PW




Foetal face



Uterus



Foetal heart




#### **UNPARALLELED EASE OF USE**

The H60 ensures ease of use with unique features such as XI STIC, Volume NT & IT $^{\text{TM}}$  and e-Motion Marker $^{\text{TM}}$ . Intuitiveness and usability are the two intelligent traits of the revolutionary H60.

#### XI STIC

XI STIC provides 3D rendered view of the foetal heart. It allows you to promptly check foetal heart beating with realistic images.



XI STIC

#### Volume NT & IT™

User-friendly Volume NT & IT™ allows improved mid-sagittal views and easier measurements. Stored volume data makes reviews and reassessments simpler.



NT measurements with Volume NT & IT™

#### e-Motion Marker™

e-Motion Marker™ displays the direction of the Transducer's beam plane on the screen and provides reference information to assist in diagnosis. The position of the uterus and ovaries can be expressed easily and intuitively.



e-Motion Marker™


 $<sup>\</sup>ensuremath{^*}\xspace$  Above features may not be available in some countries.

### **ADVANCED TECHNOLOGY THAT SAVES TIME**

QuickScan<sup>™</sup>, digital TGC preset and intelligent control panel provide the H60 with time saving functions that offer the fastest and most accurate readings possible.


#### QuickScan™

Auto-optimisation of important imaging parameters is possible with just the touch of a button, and the end result is the maximisation of workflow efficiency.









CCA Doppler with QuickScan™

#### **Intelligent Control Panel**








User customizable touch menus Scroll & click toggle button



User customisable keys



The H60's innovative features are efficiently integrated into a sleek and compact space-saving design. The built-in printer and movable base allow ease of use across multiple locations.



#### Slim & Lightweight

The slim and compact design allows the H60 to fit in a patient's bed for the ultimate convenience.

#### **Lift Control Panel**

Allows smooth movement and maximises user comfort.



## **Encased Printer & Basket Space**

The built-in printer and basket place efficiently maximise storage space.







#### **COMPREHENSIVE RANGE OF TRANSDUCERS**

Samsung Medison's comprehensive range of Transducers ensures you find the perfect device as per your patient's needs.

#### **Curved Array Transducers**



#### CA2-8AD

- Application : Abdomen, Obstetrics & Obstetrics & Gynaecology Gynaecology
- Field of View: 70°

CA1-7AD

- Application: Abdomen,
- Field of View: 58°



- Application : Abdomen, Obstetrics & Gynaecology
- Field of View: 58°



CF4-9

- Application : Paediatric & Vascular
- Field of View: 92°



- Application : Abdomen, Obstetrics & Gynaecology
- Field of View: 68°

#### **Endo-Cavity Transducers**



VR5-9

- Application : Obstetrics, Gynaecology & Urology
- Field of View: 150°

EVN4-9

- Application: Obstetrics, Gynaecology & Urology
- Field of View: 148°
- ER4-9
- Application : Obstetrics, Gynaecology & Urology
- Field of View: 148°

#### **Linear Array Transducers**



LF5-13

- Application : Small Parts, Vascular & Musculoskeletal
- Field of View: 40mm



LA3-14AD

- Application : Small Parts, Vascular & Musculoskeletal
- Field of View: 50mm





L5-13

- Application : Small Parts, Vascular & Musculoskeletal
- Field of View: 40mm

#### **Phased Array Transducers**





- Application : Cardiac, Abdomen & TCD
- Field of View: 90°

PE2-4



#### **SP3-8**

- Application : Cardiac, Abdomen & TCD
  - Field of View: 90°

#### **Volume Transducers**



CV1-8AD

• Application : Abdomen, Obstetrics & Gynaecology • Field of View: 70°

#### 3D2-6

• Application : Abdomen, Obstetrics & Gynaecology • Field of View: 84°



Obstetrics & Gynaecology • Field of View: 69.78°



#### V5-9

- Application : Abdomen, Application: Obstetrics, Gynaecology & Urology
  - Field of View: 150°



#### 3D4-9

- Application: Obstetrics, Gynaecology & Urology
- Field of View: 145°

#### **CW Pencil Type Transducers**



#### DP2B

- Application : Cardiac
- CW4.0
- Application : Cardiac

Samsung Ultrasound H60 | 15 Performance in Style